



Diamant-Scheiben galvanisch gebunden  
Diamond Discs galvanic bound



914 einseitig  
914 single-sided



915 doppelseitig  
915 double-sided



916 einseitig  
916 single-sided



916 einseitig  
916 single-sided



917 doppelseitig  
917 double-sided





918 doppelseitig  
918 double-sided

**Beschreibung:** Diamant-Scheiben galvanisch gebunden.  
**Anwendungsbereich:** Zum Trennen, Separieren und seitlichen Beschleifen von Keramik, Kunststoffen, Composite, Gips, und Metallen.

**Description:** Diamond Discs galvanic bound.  
**Field of application:** For cutting, separating, and lateral grinding of ceramic, plastic, composite, plaster and metals.

Bestellnummer / Order Number			Ø 19 mm	Ø 22 mm
Montiert mounted	unmontiert unmounted			
HP 914/190 F	914/190 F	(fein, einseitiger Belag; fine, single-sided coating)	<input type="checkbox"/>	
HP 914/220 F	914/220 F	(fein, einseitiger Belag; fine, single-sided coating)		<input type="checkbox"/>
HP 914/190 XF	914/190 XF	(x-fein, einseitiger Belag; x-fine, single-sided coating)	<input type="checkbox"/>	
HP 914/220 XF	914/220 XF	(x-fein, einseitiger Belag; x-fine, single-sided coating)		<input type="checkbox"/>
HP 915/190 F	915/190 F	(fein, doppelseitiger Belag; fine, double-sided coating)	<input type="checkbox"/>	
H 915/220 F	915/220 F	(fein, doppelseitiger Belag; fine, double-sided coating)		<input type="checkbox"/>
HP 915/190 XF	915/190 XF	(x-fein, doppelseitiger Belag; x-fine, double-sided coating)	<input type="checkbox"/>	
HP 915/220 XF	915/220 XF	(x-fein, doppelseitiger Belag; x-fine, double-sided coating)		<input type="checkbox"/>
HP 916/190 F	916/190 F	(fein, einseitiger Belag; fine, single-sided coating)	<input type="checkbox"/>	
HP 916/220 F	916/220 F	(fein, einseitiger Belag; fine, single-sided coating)		<input type="checkbox"/>
HP 916/190 M	916/190 M	(mittel, einseitiger Belag; medium, single-sided coating)	<input type="checkbox"/>	
HP 916/220 M	916/220 M	(mittel, einseitiger Belag; medium, single-sided coating)		<input type="checkbox"/>
HP 917/190 F	917/190 F	(fein, doppelseitiger Belag; fine, double-sided coating)	<input type="checkbox"/>	
HP 917/220 F	917/220 F	(fein, doppelseitiger Belag; fine, double-sided coating)		<input type="checkbox"/>
HP 917/190 M	917/190 M	(mittel, doppelseitiger Belag; medium, double-sided coating)	<input type="checkbox"/>	
HP 917/220 M	917/220 M	(mittel, doppelseitiger Belag; medium, double-sided coating)		<input type="checkbox"/>
HP 918/220 F	918/220 F	(fein, doppelseitiger Belag; fine, double-sided coating)		<input type="checkbox"/>





Sinterdiamant-Scheiben / Nonflex  
Sinterdiamond Discs / Nonflex



**Beschreibung:** Sinterdiamant-Scheiben / Nonflex

**Description:** Sinterdiamond discs / Nonflex

Bestellnummer Order Number	Durchmesser Ø Diameter Ø	Stärke Strength	Anwendungsbereich Field of application
M 941/200 HP (mittel, medium)	Ø 20 mm	0,27 mm	Keramik; ceramic
M 941G /220 HP (grob, coarse)	Ø 22 mm	0,30 mm	Metall, Keramik; metal, ceramic
M 941XG /220 HP (x-grob, x-coarse)	Ø 22 mm	0,30 mm	Metall, Keramik; metal, ceramic

Sinterdiamant-Scheiben / Mediumflex  
Sinterdiamond Discs / Mediumflex




**Beschreibung:** Sinterdiamant-Scheiben / Mediumflex

**Description:** Sinterdiamond discs / Mediumflex

Bestellnummer Order Number	Durchmesser Ø Diameter Ø	Stärke Strength	Anwendungsbereich Field of application
M 942a /210 HP (mittel, medium)	Ø 21 mm	0,20 mm	Kunststoff, Gold; plastic, gold
M 942b /210 HP (mittel, medium)	Ø 21 mm	0,23 mm	Kunststoff, Gold; plastic, gold
M 942 /300 HP (mittel, medium)	Ø 30 mm	0,30 mm	Gips; plaster
M 942 /400 HP (mittel, medium)	Ø 40 mm	0,30 mm	Gips; plaster
M 942 /450 HP (mittel, medium)	Ø 45 mm	0,30 mm	Gips; plaster

Sinterdiamant-Scheiben / Flex  
Sinterdiamond-Discs / Flex



**Beschreibung:** Sinterdiamant-Scheiben / Flex

**Description:** Sinterdiamond discs / Flex


Bestellnummer Order Number	Durchmesser Ø Diameter Ø	Stärke Strength	Anwendungsbereich Field of application
M 943F /210 HP (fein, fine)	Ø 21 mm	0,16 mm	Keramik; ceramic
M 943F /250 HP (fein, fine)	Ø 25 mm	0,16 mm	Keramik; ceramic



Sinterdiamant-Scheiben / Superflex  
Sinterdiamond-Discs / Superflex

**Beschreibung:** Sinterdiamant-Scheiben / Superflex

**Description:** Sinterdiamond discs / Superflex



Bestellnummer Order Number	Durchmesser Ø Diameter Ø	Stärke Strength	Anwendungsbereich Field of application
M 944F /080 HP (fein, fine)	Ø 0,8 mm	0,15 mm	Keramik; ceramic
M 944F /120 HP (fein, fine)	Ø 12 mm	0,15 mm	Keramik; ceramic
M 944F /160 HP (fein, fine)	Ø 16 mm	0,15 mm	Keramik; ceramic
M 944F /180 HP (fein, fine)	Ø 18 mm	0,15 mm	Keramik; ceramic
M 944F /210 HP (fein, fine)	Ø 21 mm	0,15 mm	Keramik; ceramic
M 944a /210 HP (mittel, medium)	Ø 21 mm	0,17 mm	Keramik; ceramic
M 944b /210 HP (mittel, medium)	Ø 21 mm	0,23 mm	Keramik; ceramic

Notiz / Notice: