

# Stiff Reamers



	Gesamtlänge <i>Total Length</i>	Spiral-Länge <i>Twist Length</i>	AO Art. No.	HARRIS Art. No.	HUDSON Art. No.	STRYKER / ZIMMER Art. No.
--	------------------------------------	-------------------------------------	----------------	--------------------	--------------------	------------------------------



Ø 8,0 mm	350 mm	185,0 mm	20350/080	21350/080	22350/080	23350/080
Ø 8,5 mm	350 mm	140,0 mm	20350/085	21350/085	22350/085	23350/085
Ø 9,0 mm	350 mm	158,5 mm	20350/090	21350/090	22350/090	23350/090
Ø 9,5 mm	350 mm	158,5 mm	20350/095	21350/095	22350/095	23350/095
Ø 10,0 mm	350 mm	158,5 mm	20350/100	21350/100	22350/100	23350/100
Ø 10,5 mm	350 mm	158,5 mm	20350/105	21350/105	22350/105	23350/105
Ø 11,0 mm	350 mm	158,5 mm	20350/110	21350/110	22350/110	23350/110
Ø 11,5 mm	350 mm	158,5 mm	20350/115	21350/115	22350/115	23350/115
Ø 12,0 mm	350 mm	168,5 mm	20350/120	21350/120	22350/120	23350/120
Ø 12,5 mm	350 mm	168,5 mm	20350/125	21350/125	22350/125	23350/125



Andere Größen und Anschlüsse auf Anfrage / other sizes and couplings upon request



# Stiff Reamers

Markraumbohrer

Markraumbohrer / Stiff Reamers

	Gesamtlänge <i>Total Length</i>	Spiral-Länge <i>Twist Length</i>	AO Art. No.	HARRIS Art. No.	HUDSON Art. No.	STRYKER / ZIMMER Art. No.
--	------------------------------------	-------------------------------------	----------------	--------------------	--------------------	------------------------------




Ø 13,0 mm	350 mm	168,5 mm	20350/130	21350/130	22350/130	23350/130
Ø 13,5 mm	350 mm	168,5 mm	20350/135	21350/135	22350/135	23350/135
Ø 14,0 mm	350 mm	168,5 mm	20350/140	21350/140	22350/140	23350/140
Ø 14,5 mm	350 mm	168,5 mm	20350/145	21350/145	22350/145	23350/145
Ø 15,0 mm	350 mm	168,5 mm	20350/150	21350/150	22350/150	23350/150
Ø 15,5 mm	350 mm	168,5 mm	20350/155	21350/155	22350/155	23350/155
Ø 16,0 mm	350 mm	168,5 mm	20350/160	21350/160	22350/160	23350/160
Ø 16,5 mm	350 mm	168,5 mm	20350/165	21350/165	22350/165	23350/165
Ø 17,0 mm	350 mm	168,5 mm	20350/170	21350/170	22350/170	23350/170
Ø 17,5 mm	350 mm	168,5 mm	20350/175	21350/175	22350/175	23350/175



Andere Größen und Anschlüsse auf Anfrage / other sizes and couplings upon request

# Stiff Reamers



	Gesamtlänge <i>Total Length</i>	Spiral-Länge <i>Twist Length</i>	AO Art. No.	HARRIS Art. No.	HUDSON Art. No.	STRYKER / ZIMMER Art. No.
						
Ø 18,0 mm	350 mm	168,5 mm	20350/180	21350/180	22350/180	23350/180
Ø 18,5 mm	350 mm	168,5 mm	20350/185	21350/185	22350/185	23350/185
Ø 19,0 mm	350 mm	178,5 mm	20350/190	21350/190	22350/190	23350/190
Ø 19,5 mm	350 mm	178,5 mm	20350/195	21350/195	22350/195	23350/195
Ø 20,0 mm	350 mm	178,5 mm	20350/200	21350/200	22350/200	23350/200
Ø 20,5 mm	350 mm	178,5 mm	20350/205	21350/205	22350/205	23350/205
Ø 21,0 mm	350 mm	178,5 mm	20350/210	21350/210	22350/210	23350/210
Ø 21,5 mm	350 mm	178,5 mm	20350/215	21350/215	22350/215	23350/215
Ø 22,0 mm	350 mm	178,5 mm	20350/220	21350/220	22350/220	23350/220



Andere Größen und Anschlüsse auf Anfrage / other sizes and couplings upon request